



Biochemical Test kits LabAssayTM Series

The LabAssay[™] Series is a collection of biochemical assay kits for analyzing samples from humans, mice, and other animals. Since analyses are performed in microtiter plates, many samples can be concurrently evaluated, and only small sample amounts are needed. ***For Research Use Only, Not for Diagnostic Use**

Code No.	Measuring factor	Grade	Animal species	Sample	Standard curve range	Sample amount	Measuring time	Package size
295-78901	Ammonia nitrogen		Mouse, Rat, Human	Blood	100∼400µg/dL	70µL	approx. 70min	700 tests
291-58601	ALP		Mouse, Human	Serum	0.0625~0.5mmol/L	20µL	approx. 20min	900tests
294-65801	Cholesterol		Mouse, Human	Serum	50~592.2mg/dL	2µL	approx. 10min	1,000tests
290-65901	Creatinine		Mouse, Human	Serum	2.5~10mg/dL	50µL	approx. 40min	500test
298-65701	Glucose	for Cellbilogy	Mouse, Human	Serum	50~500mg/dL	2µL	approx. 10min	1,000tests
294-63601	Free Fatty Acid (NEFA)		Mouse, Human	Serum	0.5~1.97mEq/L ※Oleic acid 1mEq=1mmol	4µL	approx. 20min	750tests
296-63801	Phospholipid		Mouse, Human	Serum	150~596.1mg/dL	2µL	approx. 10min	1,300tests
290-63701	Triglycelide		Mouse, Human	Serum	100~888mg/dL	2µL	approx. 10min	1,000tests

Note; 294-65801 LabAssay[™] Cholesterol is not available for sale in the US. 290-65901 LabAssay[™] Creatinine is not available for sale in Europe. 294-63601 LabAssay[™] NEFA is not available for sale in the US and Europe. 296-63701 LabAssay[™] Phospholipid is not available for sale in the US.

LabAssay[™]Ammonia

Ammonia is mainly produced in the intestine or the kidney. Ammonia produced in the body is converted to urea through a series of reactions known as the urea cycle in the liver and eliminated through urine.

Assay Principle (Modified Fujii-Okuda Method)

Ammonia is converted to Dioxydiphenylamine by the reaction of Phenol and Sodium Pentacyanonitrosyl ferrate(III). A reaction of the Dioxydiphenylamine with Sodium Hypochlorite produces Indophenol ,which pigments blue, under alkaline conditions. LabAssay ™ Ammonia is a kit used for the quantitative determination of ammonia nitrogen in samples by measuring absorbance of the blue color.

Kit Contents

1	Deproteinizing	100mL×1vial
2	Chromogen Reagent A	50mL×1vial
3	Chromogen Reagent B	25mL×1vial
4	Chromogen Reagent C	50mL×1vial
5	Ammonia Standard Solution	15mL×1vial
6	Dilute olution for Standard	20mL×1vial

Performance

Standard curve range
Measurement time

e 100~400µg/dL(µg/100mL) approx. 70 minutes

70µL

- Amount of sample
- Measurement wavelength 630nm





Code No.	Product name	Grade	Package size	Storage Condition
295-78901	LabAssay [™] Ammonia	for Cellbiology	700 tests	Keep at 2-10°C

LabAssay[™]ALP

Alkaline Phosphatase (ALP) is an enzyme distributed in a variety of tissues such as the liver, bone and small intestines in animals. Especially, it is used one of the Osteogenesis Markers in bone metabolism research area.

Assay principle (Alkaline Phosphatase activity assay with p-Nitrophenyl Phosphate as a substrate)

p-Nitrophenylphosphate is hydrolyzed into p-Nitrophenol in the presence of Alkaline Phosphatase(ALP). LabAssay ™ ALP is a kit designed to determine Alkaline Phosphatase activity in samples by measuring the amount of p-Nitrophenol released by p-Nitrophenylphosphate as a phosphatase substrate.

Kit Contents

1	Substrate Tablet	20 tablets
2	Buffer Solution	100mL×1vial
3	Stop Solution	100mL×1vial
4	Standard Solution	10mL×1vial

Performance

•	Standard curve range	0.0625~0.5mmol/L
•	Measurement time	approx. 20 minutes
•	Amount of sample	20µL

Measurement wavelength 405nm

Standard Curve



[Reference]

1) Ito, S. et al. : J. Pharmacol. Exp. Ther., 333, 341(2010). XExtraction liqiud of mouse kidney 2) Matsuyama, A. et al. : Clin. Exp. Pharmacol. Physiol., 45, 75(2018). XMC3T3-E1cell, C2C12Ccell 3) Chiba, T. et al. : J. Atheroscler. Thromb, 23, 1099(2016). XMouse plasma

4) Kohno, Y. et al. : Stem Cell Res. Ther., 8, 115(2017). *Mesenchymal stem cell

Code No.	Product name	Grade	Package size	Storage Condition
291-58601	LabAssay [™] ALP	for Cellbiology	900 tests	Keep at 2-10℃

LabAssav[™]Cholesterol

Cholesterol, a major component of cell membranes and the starting material in steroid synthesis in many animals, is a factor behind arteriosclerosis and other vascular diseases.

Assay Principle (Cholesterol Oxidase · DAOS method)

Hydrogen peroxide produced by a reaction of cholesterol and cholesterol oxidase let N-Ethyl-N- (2-hydroxy-3sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) and 4-Aminoantipyrin oxidize and condensate. LabAssay ™ Cholesterol is a kit to determine total cholesterol by measuring absorbance of the blue color which is generated by the oxidative condensation reaction.

Kit Content

1	Buffer	150mL×2vials
2	Chromogen	for 150mL×2vials
3	Standard Solution	10mL×1vial

Performance

Standard curve range

· Measurement time

- 50~592.2mg/dL(mg/100mL) approx. 10 minutes 2µL
- · Amount of sample Measurement wavelength



Standard Curve



[Reference]

2) Gao, F. et al. : Evid. Based Complement. Alternat. Med., 2015, 801291(2015). * Extraction liquid of rat kidney

3) Yoshioka, H. and Onosaka, S. : Fundam. Toxicol. Sci., 3, 151(2016). **Mouse plasma

4) Fujii, N. et al. : Aging Cell, 16, 508(2017). XRat plasma

Code No.	Product name	Grade	Package size	Storage Condition	
294-65801	LabAssay [™] Cholesterol	for Cellbiology	1,000 tests	Keep at 2-10℃	
*Not available for sale in the US.					

LabAssay[™]Creatineine

Creatinine is a metabolite which is produced by creatine phosphate directly or dehydration of creatine in muscle and nerve. It is excreted from the body through kidney glomerular filtration.

Assay Principle (Jaffé method)

LabAssay ™ Creatinine can be used to measure the creatinin levels in samples by Jaffe's reaction where creatinine produces quantitatively an orange color with picric acid in alkaline condition.



[Reference]

Code No.	Product name	Grade	Package size	Storage Condition
290-65901	LabAssay [™] Creatinine	for Cellbiology	500 tests	Keep at 2-10℃

XNot available for sale in Europe.

LabAssay[™]Glucose

Sugar is one of the most important sources of energy in biology. It is regulated by various factors within an organism. Glucose converges to a stable ratio of α -form and β -form in solutions.

Assay Principle (Mutarotase · GOD method)

α-D-Glucose is converted to β-D-Glucose by mutarotase. Hydrogen peroxide, which is produced by a reaction between β-D-Glucose and glucose oxidase(GOD), promotes oxidative condensation of phenol with 4-aminoantipyrine quantitatively. LabAssay TM Glucose is a kit used for the quantitative determination of glucose concentrations in samples by measuring absorbance of a red color which is generated by the oxidative condensation reaction.

Kit Contents

1	Buffer	150mL×2vial
2	Chromogen Reagent	for 150mL×2vials
3	Glucose Standard I	10mL×1vial
4	Glucose Standard II	10mL×1vial

Performance

- Standard curve range
 - Measurement time approx
- Amount of sample

2.5~10mg/dL(mg/100mL) approx. 40 minutes 50µL

Measurement wavelength 505nm(Main), 600(Sub)

Standard Curve



[Reference]

1) Yamashita, Y. et al.: Biosci. Biotechnol. Biochem., 77, 888 (2013). **Mouse plasma

3) Yamasaki, M. et al. : Food Sci. Technol. Res., 21, 827(2015). **Mouse serum

Code No.	Product name	Grade	Package size	Storage Condition
298-65701	LabAssay [™] Glucose	for Cellbiology	1,000 tests	Keep at 2-10℃

LabAssay[™]NEFA

NEFA (Non - esterified fatty acid) in the blood is transported complexed with an albumin to peripheral tissues. They are important sources of fuel for the peripheral tissues. The concentration of NEFA in the blood is regulated by a release from the adipose tissues, a consumption in the periphral tissues or a take up from the liver.

Assay Principle (ACS · ACOD method)

NEFA (Non - esterified fatty acid) forms Acyl-CoA in the presence of Acyl-CoA synthetase(ACS). Hydrogen peroxide. which is produced by a reaction between the Acyl-CoA and Acyl-CoA oxidase(ACOD), promotes oxidative condensation of 3-methyl-N-ethyl-N (-β-hydroxyethyl) -aniline (MEHA) with 4-aminoantipyrine. LabAssay ™ NEFA can be used for the quantitative determination of NEFA in samples by measuring absorbance of a blue purple color which is generated by the oxidative condensation reaction.

Kit Contens

1	Chromogen Reagent A	for 10mL×6vials
2	Solvent A	65mL×1vial
3	Chromogen Reagent B	for 20mL×6vials
4	Solvent B	130mL×1vial
5	NEFA Standard Solution (Oleic acid 1mEq/L)	10mL×1vial
Perfo	XOleic :	acid 1mEq=1mmol
Standard curve range		~1.97mEa/L



Amount of sample Measurement wavelength

Measurement time

550nm

[Reference]

2) Gao, F. et al. : Evid. Based Complement. Alternat. Med., 2015, 801291(2015). * Extraction liquid of rat kidney

approx. 20 minutes

4µL

3) Ogawa, K. et al. : Reprod. Med. Biol., online(2018). doi.org/10.1002/rmb2.12084 % Follicular fluid derived from pig vesicular ovarian follicle 4) Wang, F. et al.: J. Mol. Endocrinol., 52, 133 (2014). ** Mouse plasma

Code No.	Product name	roduct name Grade		Storage Condition	
294-63601	LabAssay [™] NEFA	for Cellbiology	750 tests	Keep at 2-10℃	

XNot available for sale in the US and Europe.

LabAssay[™]Phospholipid

Phospholipids are known as not only as major component of cell membranes but also perform vital functions within the body such as emulsification and absorption of fats or coagulation of blood.

Assay Principle (Choline Oxidase · DAOS method)

Phospholipids are hydrolyzed by Phospholipase D to release hydrogen peroxide. The formed hydrogen peroxide promotes oxidative condensation of N-ethyl-N- (2-hydroxy-3-sulfopropyl) -3,5-dimethoxyaniline sodium salt (DAOS) with 4aminoantipyrine. LabAssay ™ Phospholipid is a kit to determine Phospholipid concentration in samples by measuring absorbance of a blue color which is generated by the oxidative condensation reaction.

Kit Contens

1	Buffer	50mL×8vials
2	Chromogen Substrate	for 50mL×8vials
3	Standard Solution	10mL×2vials

Performance

- Standard curve range
- Measurement time Amount of sample

150~596.1mg/dL(mg/100mL) approx. 10 minutes 2µL 600nm(Main), 700nm(Sub)

Measurement wavelength

Standard Curve



[Reference]

- 1) Tatematsu, Y. et al. : Biol. Pharma. Bull., 41, 319(2018). XLiposome 2) Kuge, H. et al.: J. Biol. Chem., 289, 26783 (2014). XLiposome
- 3) Kessler, E. C. et al.: J. Dairy. Sci., 97, 5481 (2014). **Bovine

4) Xu, Q. et al.: Biosci. Biotechnol, Biochem., 77, 1390 (2013). * Extraction liquid of mouse kidney

Code No.	Product name	Grade	Package size	Storage Condition
296-63801	LabAssay [™] Phospholipid	for Cellbiology	1,300 tests	Keep at 2-10°C

Not available for sale in the US.



Standard Curve

LabAssay[™]Tryglyceride

Triglycerides are neutral fats consisting of three fatty acids esterified to a glycerol backbone. There are triglycerides, cholesterol, phospholipids, free fatty acids and fat-soluble vitamins as lipid-soluble substances in the blood.

Assay Principle

Triglycerides are converted to glycerol-3-phosphate by lipoprotein lipase and glycerolkinase. Hydrogen peroxide, which is produced by a reaction between the glycerol-3-phosphate and glycerol-3-phosphate oxidase(GPO), promotes oxidative condensation of N-ethyl-N- (2-hydroxy-3-sulfopropyl) -3,5-dimethoxyaniline sodium salt(DAOS) with 4-aminoantipyrine. LabAssay [™] Triglyceride can be used to detect triglycerides concentration in samples by measuring absorbance of a blue color which is generated by the oxidative condensation reaction.

Kit Contents

1	Buffer	105mL×3vials
2	Chromogen Substrate	for 105mL×3vials
3	Standard Solution	10mL×1vial

Performance

Standard curve rangeMeasurement time

Amount of sample

100~888mg/dL(mg/100mL) approx 10 minutes

- Measurement wavelength
- approx. 10 minutes 2µL 600nm(Main), 700nm(Sub)



[Reference]

1) Gao, F. et al. : Evid. Based Complement. Alternat. Med., 2015, 801291(2015). * Extraction liquid of rat kidney

2) Funakoshi, T. et al. : Biochem. Biophys. Rep., 13, 39(2018). XRat primary muscle satellite cell

- 3) Moser, V. A. and Pike, C. J. : eNeuro, 4, e0077-17(2017). **Mouse plasma
- 4) Fujii, N. et al. : Aging Cell, 16, 508(2017). XRat plasma

5) Oh, T-W., et al. : J. Exerc. Nutrition Biochem., 19, 247(2015). XRat plasma

6) Fan, Y. et al. : J. Biomed. Sci., 23, 56(2016). * Mouse serum

Code No.	Product name	Grade	Package size	Storage Condition
290-63701	LabAssay [™] Triglyceride	for Cellbiology	1,000 tests	Keep at 2-10℃

Reagents for Adiposity & Diabetes Research

Wako offers a wide variety of products for studying Obesity and Diabetes.



Code No.	Product name	Standard curve range	Grade	Package Size	Storage condition
299-75501	GLP-1 ELISA Kit Wako , High Sensitive	0.94~30pmol/L	for Diabetes Research	96 tests	Keep at 2-10℃
		(3.1~100pg/mL)	IOI Diabetes Research	50 (0303	
202 70201	Active GLP-1 ELISA Kit Wako,	0.12330pmol/l	for Diabotos Posoarch	OC to ata	Keep at 2-10℃
293-79301	Chemiluminescent	0.123/03001101/2	IOI Diabetes Research	90 lesis	
291-59201	Rat GLP-1 ELISA Kit Wako	206~50,000pg/mL	for Diabetes Research	96 tests	Keep at 2-10℃
292-60601	Rat GLP-2 ELISA Kit Wako	0.137~100ng/mL	for Diabetes Research	96 tests	Keep at 2-10℃
299-73801	Mouse GIP(Active) ELISA Kit Wako	7.8~500pg/mL	for Diabetes Research	96 tests	Keep at 2-10℃
292-80001	Glucagon ELISA Kit Wako	2.2~143.6pmol/L(7.8~500pg/mL)	for Diabetes Research	96 tests	Keep at 2-10℃
291-73501	Mouse/Rat PYY ELISA Kit Wako	0.15~12.5ng/mL	for Adiposity Research	96 tests	Keep at 2-10℃
295-57401	Rat C-Peptide ELISA Kit Wako	1.56~50ng/mL	for Diabetes Research	96 tests	Keep at 2-10℃
297-57601	Rat Leptin ELISA Kit <i>Wako</i>	78.1~5,000pg/mL(Non-Plasma、Non-serum)	for Adiposity Research	96 tests	Keep at 2-10℃
		312.5~20.000pg/mL(Plasma, Serum)			

LBIS

ELISA for Measuring Insulin

Code No.	Maker Code	Product name	Standard curve range	Package Size	Price
637-10629	AKRIN-010H	Lbis Insulin-Rat (H type)	$0.5{\sim}100$ ng/mL	96 tests	Keep at 2-10℃
634-10379	AKRIN-011H	Lbis Insulin-Mouse (H type)	$0.5{\sim}100$ ng/mL	96 tests	Keep at 2-10℃
631-01479	AKRIN-010T	Lbis Insulin-Rat-T	$0.156{\sim}10$ ng/mL	96 tests	Keep at 2-10℃
638-01489	AKRIN-011T	Lbis Insulin-Mouse-T	$0.156{\sim}10$ ng/mL	96 tests	Keep at 2-10℃
631-07199	AKRIN-010S	Lbis Insulin-Rat (S type)	$0.1 \sim \! 10$ ng/mL	96 tests	Keep at 2-10℃
630-07289	AKRIN-011S	Lbis Insulin - Mouse (S type)	78 ~5,000pg/mL	96 tests	Keep at 2-10℃
Inquiry	AKRIN-130	Lbis Insulin-Rat (U-E type)	39 ~2,500pg/mL	96 tests	Keep at 2-10℃
Inquiry	AKRIN-031	Lbis Insulin-Mouse (U type)	39 ~2,500pg/mL	96 tests	Keep at 2-10℃
630-24149	AKRIN-010RU	Lbis Rat Insulin ELISA KIT(RTU)	100~12,000pg/mL	96 tests	Keep at 2-10℃
633-23919	AKRIN-011RU	Lbis Mouse Insulin ELISA KIT(RTU)	100~12,000pg/mL	96 tests	Keep at 2-10℃
637-01459	AKRIN-012T	Lbis Insulin-Dog-T	0.188~12ng/mL	96 tests	Keep at 2-10℃
634-01469	AKRIN-013T	Lbis Insulin-Porcine	$0.188{\sim}12$ ng/mL	96 tests	Keep at 2-10℃
Inquiry	AKRIN-014T	Lbis Insulin-Monkey-T	$0.156{\sim}10$ ng/mL	96 tests	Keep at 2-10℃

Other products : Reagent for Adiposity & Diabetes Research

Code No.	Maker Code	Product name	Standard curve range	Package Size	Price
630-23049 AKMPI-111		Lbis Proinslin Mouse/Rat	1.47~94.3pg/mL(Standard curve range)	60 tosts	Keep at 2-10℃
			7.35~471.5pg/mL(Amount of a sample:10µL)	ou lesis	
625-07220		Lbis C-Peptide-Mouse (U type)	46.9~3,000pg/mL(Standard curve range)	06 tosts	Keep at 2-10℃
033-07239	AKKCF-031		234.5~15,000pg/mL(Amount of a sample:10µL)	90 lesis	
622-07270		Lbis C-Peptide Rat (U type)	46.9~3,000pg/mL(Standard curve range)	06 tosts	Keep at 2-10℃
033-07279	AKKCF-030		234.5~15,000pg/mL(Amount of a sample:10µL)	90 lesis	
621-10290		Lbis Leptin-Mouse	103~25,000pg/mL(Standard curve range)	06 tosts	Keep at 2-10℃
051-10589	AKKLF-011		20.6~5,000pg/mL(Amount of a sample:10µL)	90 lesis	
638-04309	AKRAL-121	Lbis Albumin Mouse ELISA Kit	50~1,000ng/mL(Standard curve range)	96 tests	Keep at 2-10℃
635-04319	AKRAL-120	Lbis Albumin Rat ELISA Kit	50~1,000ng/mL(Standard curve range)	96 tests	Keep at 2-10℃
			31.3~2,000pg/mL(Standard curve range)		
639-13749	AKRGH-010	Lbis GH-Rat	313~20,000pg/mL(Amount of a sample:5µL)	96 tests	Keep at 2-10℃
			62.6~4,000pg/mL(Amount of a sample:25µL)		
			313~10,000pg/mL(Standard curve range)	06 tosts	Koop at 2, 10°C
630-23929	AKKLH-0105	LDIS KALLH ELISA KIT(S type)	1.56~50ng/mL(Amount of a sample:10µL)	96 lests	Keep at 2-10°C
(22, 21720		AKRTS-010S2 Lbis TSH-Rat (S II type)	184~18,000pg/mL(Standard curve range)	OC hosts	Kana at 2, 10%
033-31729	AKK15-01052		0.92~90ng/mL(Amount of a sample:20µL)	90 lesis	Keep at 2-10 C
			1.56~50pg/mL(Standard curve range)		
637-15129	AKMGP-011	Lbis GLP-1(active) ELISA KIT	7.8~250pg/mL(Amount of a sample:10µL)	96 tests	Keep at 2-10℃
			3.9~125pg/mL(Amount of a sample:20µL)		
638-13079		KMAN-011 Lbis High Molecular Adiponectin-Mouse/Rat	3.13~200ng/mL(Standard curve range)	96 tests	Keep at 2-10℃
030-130/9			78.25~5,000ng/mL(Using sample diluted 25times)	30 10515	

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from Wako; http://www.e-reagent.com / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Pure Chemical Corporation

http://ffwk.fujifilm.co.jp/en/index.html 1-2, Doshomachi 3-Chome Chuo-Ku, Osaka 540-8605, Japan Tel.: +81 6 6203 3741 Fax: +81 6 6203 1999 ffwk-cservice@fujifilm.com Online Catalog: www.e-reagent.com

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road,Richmond, VA 23237, U.S.A. Toll-Free (U.S. only): +1 877 714 1920 Tel.: +1 804 714 1920 Fax: +1 804 271 7791 wkuslabchem@fujiflim.com

FUJIFILM Wako Chemicals Europe GmbH

www.wako-chemicals.de Fuggerstr. 12 D-41468 Neuss, Germany Tel.: +49 2131 311 0 Fax: +49 2131 311 100 labchem@wako-chemicals.de